





AI & EEG

In our project, we used Machine Learning, to enhance EEG analysis by:

- •Learning complex patterns from preprocessed EEG data that are difficult to detect manually.
- •Classifying brain states such as normal and abnormal activity based on EEG signal features.
- •Improving diagnostic assistance, particularly in tasks like:
 - Abnormal EEG detection
 - Supporting neurological screening
 - Enabling future integration with assistive technologies (e.g., brain-controlled systems)

System Architecture

Frontend:

Built with HTML/CSS and JavaScript frameworks for interactivity.

Backend:

Developed using Flask (Python) to handle server-side operations and cloudflare database.

Machine Learning Model:

Utilizes a Multi-Layer Perceptron (MLP) neural network trained to classify EEG signals.

Data Handling:

EEG data is preprocessed and scaled before being fed into the model for prediction.

Machine Learning Methodology

1. Dataset & Labeling:

- •We utilized the **Temple University Hospital Abnormal EEG Corpus (TUAB v3.0.1)** a clinically validated dataset designed for EEG signal classification.
- •Access was obtained by submitting a **Data Use Agreement (DUA)** through the official TUH EEG portal

Dataset Composition

- •EEG recordings are stored in **European Data Format (EDF)**.
- •Each file contains:
 - Multichannel brainwave signals from up to 36 scalp electrodes
 - Patient metadata (age, gender, recording date)
 - Sampling frequencies typically between 250–500 Hz

MNE?

MNE is an open-source Python library specialized for processing, analyzing, and visualizing **EEG** data.

It supports reading **EDF files**, applying filters, handling channels, and extracting features for machine learning.

EEG signals are loaded using the MNE Python library, which provides tools such as read_raw_edf() for reading files, filter() for bandpass filtering, and get_data() for accessing signal arrays — enabling efficient preprocessing and analysis.

- •mne.io.read_raw_edf() to read EDF EEG files
- •raw.filter(l_freq, h_freq) to apply bandpass filter
- raw.get_data() to extract signal data as NumPy arrays

Data Representation

- EEG signals are loaded using the MNE Python library, producing:
 - EEG data is stored as a 2D floating-point array with shape (channels × samples).
 For example, (36, 10,000) means the recording has 36 EEG channels, each containing 10,000 time points.

The number of samples represents how many times the signal was recorded per second — this is controlled by the **sampling rate**.

If the sampling rate is **250 Hz**, it means the signal was recorded **250 times per second**, so **10,000** samples ≈ **40 seconds** of data.

This format allows us to apply filters, extract features, and segment the data for machine learning.

- Each channel is a 1D NumPy array of type float64, representing voltage values in microvolts (μV)
- These continuous signals are segmented into fixed-length windows (e.g., 1000 samples) for model input.

Data preprocessing

- Each EEG file is labeled as:
 - 0 → Normal
 - 1 → Abnormal
 - based on folder path and expert annotations.
- Every extracted window inherits the label of its parent recording.

Artifact Removal

- EEG recordings often include **non-brain signals** known as **artifacts**, which interfere with true brainwave patterns.
- Common sources of artifacts:
 - Muscle movement (e.g., jaw clenching)
 - Eye blinks and body motion
 - Loose or poorly connected electrodes
- Since these artifacts can distort the signal and confuse the model, we applied **automatic artifact detection** only on **normal EEG recordings** to avoid accidentally removing pathological information.
- We used the annotate_muscle_zscore() function from the MNE library, which:
 - Analyzes the EEG signal in the 30–90 Hz range
 - Tags segments with strong muscle-related noise
 - Automatically removes those segments from the data

Filtering and Scaling

- Following artifact removal, we performed frequency filtering to isolate brainwave ranges relevant for analysis.
- A bandpass filter (0.5–40 Hz) was applied to:
 - Eliminate low-frequency drifts (<0.5 Hz), often caused by motion or sweat
 - Remove high-frequency noise (>40 Hz), including electrical interference and residual muscle signals
 - Retain only clinically meaningful frequencies (delta, theta, alpha, beta)
- Unlike artifact removal, which removes specific segments, filtering adjusts the full signal to clean its frequency content.
- We also converted EEG signals from Volts to Microvolts (μ V) to meet clinical standards using a scaling factor of 1e6.
- These operations were performed using the Preprocessor() class from Braindecode, which:
 - Applies bandpass filtering
 - Performs unit conversion
 - Ensures uniform preprocessing across all recordings
- The result is a set of standardized, noise-reduced EEG signals ready for segmentation and feature extraction.

Segmentation (Windowing)

- After preprocessing, each EEG recording was split into **non-overlapping windows of 1000 samples** using Braindecode's create_fixed_length_windows() function.
- This standardizes input size across all recordings, regardless of their original duration.
- Start offset was set to 0, and windows that didn't fit exactly were discarded (drop_last_window=True).
- Each window inherited the binary label (0 = normal, 1 = abnormal) from its parent EEG file.
- This step created a uniform dataset structure, making it suitable for deep learning and efficient training.

Class Balancing

- •After segmentation, the number of **abnormal windows was significantly higher** than normal ones, which could bias the model.
- •To address this, we used RandomOverSampler from imblearn to duplicate samples from the minority class (normal) and balance both classes.
- •Before applying oversampling, the window data was **flattened to 2D** (samples, features)using reshape() to match the expected input format for the sampler.
- •After resampling, the data was reshaped back into its original **3D format (samples × channels × time)** for model compatibility.

This ensured that both normal and abnormal EEG patterns were equally represented during training, improving model fairness and performance.

Handcrafted Feature Extraction

- After balancing, each EEG window (shape: channels × time) was processed to extract handcrafted features that describe the signal behavior numerically.
- •We extracted features only from those selected in the selected_feature_names.txt file to save time and reduce overfitting.
- •Features were computed for each channel independently and then concatenated into a single vector per window.

- Time-Domain Features
- → Captures basic statistical properties of the raw signal:
- Mean, Standard Deviation, Minimum, Maximum
- Skewness (asymmetry), Kurtosis (peakedness)

Feature Extraction

Frequency-Domain Features (PSD)

- •EEG signals can be analyzed by their frequency content, We used Welch's method (a fast Fourier-based approach) to calculate the average power in standard brainwave bands:
 - **Delta (0.5–4 Hz)** deep sleep
 - Theta (4–8 Hz) drowsiness or meditation
 - Alpha (8–13 Hz) calm, relaxed state
 - Beta (13–30 Hz) active thinking
 - Gamma (30–40 Hz) high-level cognitive processing

This helps identify which brain states are dominant in the signal.

Hjorth Parameters

These are three statistical measures based on signal derivatives (how fast the signal changes):

- **1.Activity** Variance of the signal → how strong or active the signal is overall
- **2.Mobility** Standard deviation of the signal's first derivative → how quickly it changes
- **3.Complexity** Compares the first and second derivatives → how "irregular" or unpredictable the signal shape is

They offer a quick snapshot of signal dynamics and structure.

Catch22 Features

- •Catch22 is a toolkit that computes 22 standardized statistical features from any time-series signal.
- •It includes:
 - Entropy how random the signal is
 - Autocorrelation how much the signal repeats over time
 - Outlier detection checks for extreme spikes
 - Nonlinearity metrics captures irregular, chaotic behavior

These features summarize a wide range of behaviors in a compact format, useful for machine learning.

Feature Selection

After extracting a large number of features from each EEG window, not all of them were useful for classification. To avoid training on irrelevant or noisy features, we applied a feature selection strategy based on model performance.

How it Worked

- •We trained an initial version of the model using all features.
- •We then measured the importance of each feature by checking how much it contributed to improving the model's F1-score (a balance between precision and recall).
- •Based on this analysis, we kept only the most impactful features and discarded the rest.

Final Setup

- •The names of the selected features were stored in a file: selected_feature_names.txt
- •During feature extraction, we computed only the selected features, skipping everything else.

 This reduced the dataset size, increased training speed, and improved model accuracy by avoiding overfitting.

Model

Train-Validation Split

- •EEG feature data and labels were split into 80% training and 20% validation sets.
- •Stratified split was used to maintain class balance.
- •Ensures that the model is trained on diverse data and evaluated fairly.

Data Loaders

- Created Data Loader objects for batching and shuffling:
 - •Batch size = 64
- •Improves training performance and generalization.

Class Weights Handling

Even after oversampling, slight class imbalance can remain due to the train-validation split. To address this, we:

- Computed final class weights using compute_class_weight() on the training set only.
- Applied these weights in CrossEntropyLoss() to ensure fair learning across both classes.
- Class weights balance the model's attention during training.

MLP Model Architecture

- •Implemented a custom Multi-Layer Perceptron (MLP) with:
 - 3 hidden layers (256 → 128 → 64 neurons)
 - BatchNorm, ReLU activations, Dropout
 - Output layer: 2 neurons (normal vs abnormal)
- Designed to balance accuracy and computational efficiency
- Optimizer: Adam with learning rate = 1e-3
- Loss Function: Weighted Cross Entropy.
- •Trained for up to 100 epochs with performance tracking.